


**#4010-02 Dual Leg Anesthesia Shield**


Two leg mounting allows unit to be placed at edge of table.



**#4010-02 Dual Leg Anesthesia**



## #4010-02 Dual Leg Anesthesia Shield Specifications



### When ordering:

Note the legs of the shield mount into rectangular bed rail holders. The legs have permanent stops to prevent slippage.

**Most** OR tables are similar in width. Please confirm dimension of table and shield width. These must be within 1/2" in order for the shield to fit properly. Pemco can manufacture the shield to the required width if needed.

